DE

EN

FR

ES

PT


IT

NO

DA

PARKINSONSTOCK

PARKINSON'S CANE

ANTI - FREEZING
FUNCTION

Art. 3200

Art. 3201

Art. 3203

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DIN EN ISO 13485
REG.-NR. MCI 0115318

This walking stick is an aid for Parkinson's patients. As Parkinson's disease progresses, sudden movement blockages can occur. The patient has the feeling of „freezing“ from one moment to the next. This phenomenon is called freezing. The laser beam is an obstacle-free, visual stimulus and serves to overcome the blockage of movement by walking towards it.

The following instructions provide you with descriptions on how to set up your Parkinson's cane correctly, explanations on how to use it properly and tips on how to enjoy your Parkinson's cane for a long time.

1. Intended use and indication	10
2. Contraindication	10
3. General safety information	10
4. Adjusting the walking aid	10
5. Use	12
6. Precautionary measure	12
7. Safe use of laser sources	13
8. Storage and care	13
9. Limitation of use	13
10. Re-use (rental) / duration of use	14
11. Warranty and liability	14
12. Disposal	14
13. Technical data	15
14. Replacement parts	16
15. Symbol directory	16

1. Intended use and indication

This walking stick is an aid for Parkinson's patients. It serves to compensate for a disability and compensates for the loss or restriction of physical functions in the execution of recurring movements and activities in everyday life. Our walking sticks are designed for people with light mobility impairment with impairment to movements, compromised balance with sufficiently intact walking ability and impaired coordination. They are used to maintain, promote and stabilise walking or to partially relieve the strain on the lower limbs. The walking sticks can be used both indoors and outdoors. **It is imperative that the safety instructions are also observed!**

2. Contraindication

The use of walking sticks is unsuitable if there is a loss of limb in both arms or joint contractures/joint damage in both arms.


3. General safety information

Do not hang any bags, pouches or other objects on the walking stick! The walking stick may under no **circumstances be used for any other purpose than that for which it is intended**, e.g. for lifting or moving objects or as a tool for switching light switches or pressing buttons on lifts! It is essential that the ground has a good grip to ensure a secure hold of the walking stick! Avoid wet floors or surfaces such as gravel, sand, foliage, ice or snow! **Pay attention to stumbling points** such as unevenness, holes in the floor, loose floor coverings or edges, because the walking stick can unexpectedly slip away or sink in! Extra caution is recommended when lighting conditions are insufficient, since irregularities in the ground can then only be seen with difficulty! When adjusting the push-button, it must protrude from the hole as far as it will go (Fig. 1)! The holes for adjusting the walking stick undergo a certain amount of wear. It is important to ensure that they are not worn out, as this would prevent the push button from locking securely. Prior to each use, make sure that all screw connections are firmly tightened (Fig. 2)!

4. Adjusting the walking stick

Warning! Please contact your specialist dealer to ensure proper functional and ergonomic adjustment! Incorrectly adjusted walking sticks may restrict the functionality of the walking stick or lead to overloading or incorrect loading of the patient! The walking stick is handed over to you as a patient in a proper condition, whereby it is individually adjusted to your body by your specialist dealer or health care specialist.


If you have any doubts as to whether this is the case, please contact the responsible specialist personnel. The walking stick can only offer you the relief you desire if it is optimally adapted to your body size.

Adjustment: The walking stick is equipped with a union nut, which must be loosened by approx. half a turn anticlockwise prior to adjustment (Fig. 2). Unlock the height adjustment by pressing in the push button. Adjust the length of the walking stick to your body size by pulling the support tubes out or pushing them together. Lock the height adjustment once again by aligning the tubes as such, that the two locking holes are on top of each other, whereby the push-button engages again. **Following adjustment, the union nut must be tightened again. Caution:** Risk of falling due to walking sticks buckling when the height adjustment is not locked correctly! Check the correct locking after each adjustment by attempting to twist or move the walking stick!

Activating the laser: The laser is activated by a button located in the module on the tube, below the grip. When the laser is switched on, the LED ring of the switch-on button lights up. Only when the walking stick is placed on the ground, the walking stick projects the laser beam onto the ground. If the walking stick is lifted, it disappears again. The laser function can be activated or deactivated at any time. This does not affect the use of the walking stick. The battery life is approx. 8 hours when the laser is in continuous operation.

Charging the walking stick: The power supply is provided by a rechargeable battery. The walking stick can be charged with the included power adapter and Magnetsafe micro USB charging cable. If the laser beam stops working, charge the walking stick. A glow of the power button is not an indication of a sufficiently charged walking stick. After approx. 2.5 hours of charging, the battery is fully charged. During the charging process, a lamp on the side of the module lights up red. When the battery is fully charged, the red lamp goes off and a blue lamp lights up.

Note: The lamp only lights up blue when the charging cable is plugged in.

5. Use

There are many different ways to use the walking sticks and it all depends on the individual restrictions of the walking ability. Seek instruction in the use of the walking stick from your specialist dealer or health care specialist! Walking sticks must not be used in water. The walking stick meets the requirements of IP protection class IP23. This means that the walking stick is protected against dust and dripping water. Permanent contact of the walking stick with moisture should nevertheless be avoided. The walking stick must be dried appropriately when it comes into contact with water. The correct handling of the walking stick takes some practice. We therefore recommend getting additional support from an assistant at the beginning. Always use the walking stick on the side opposite the injured leg or the leg that needs to be protected! Guide the walking stick close to the body, not at an angle facing outwards (Fig. 3)! You might fall otherwise! During use, grasp the handle of the walking stick completely and firmly with your hand! Always place the walking stick as straight as possible on the ground so that you do not slip away! Find chairs with armrests to sit on. You will then be able to support yourself on the armrests when standing up! Wear shoes that are as comfortable as possible with flat soles so that your leg can better withstand additional strain! **If you do fall:** If possible, let the walking stick fall to the side so that you do not fall onto it! Have someone accompany and support you if possible, so that they can catch you if you fall! Try and arrange for qualified personnel to instruct you in climbing stairs with a walking stick! When climbing stairs with a walking stick, it is important to hold on to the railing using your free hand (Fig. 4)! When climbing stairs with two walking sticks, hold both walking sticks together in one hand and make sure to hold on to the railing firmly using your free hand (Fig. 5)!

6. Precautionary measure

Using a walking stick places an unusual amount of strain on individual parts of the body! This may result in the following undesirable side effects: The hands may encounter pressure or chafing! Keep in mind that this can lead to pain and overexertion at the beginning or during prolonged exertion! It is important to take sufficient breaks!


Fig. 3


Fig. 4

Fig. 5

7. Safe use of laser sources

Do not point the laser beam at other persons! Do not look into the direct or reflected beam! If the laser beam hits the eye, deliberately close your eyes and turn away! Do not use optical instruments (e.g. magnifying glasses, binoculars) to observe the laser source! Never manipulate the laser source!

8. Storage and care

Never store the walking sticks in the vicinity of heat sources such as furnaces, heaters or in places susceptible to heat in a car! Do not expose walking sticks to sunlight for prolonged periods of time! Strong sun or heat damages the plastic! Protect your walking sticks from toppling over and from damage to the paintwork (e.g. by using straps or holders). **The rubber tip is a wearing part and must be replaced at regular intervals.** Walking sticks must only be used again after an extended period of non-use following prior inspection by qualified personnel!

Warning! Never clean the grip and rubber tip using oily substances! There is a risk of slipping! In the event of contamination, clean the walking stick using a mild cleaning agent and a soft cloth where necessary. Never use cleaners containing solvents, since these will damage the plastic! **Note:** Should disinfection be required, then a suitable disinfectant agent should be used. Please follow the manufacturer's instructions for use.

9. Limitation of use

The walking stick must only be subjected to loads in accordance with the manufacturer's instructions. **Observe the maximum permissible body weight that is specified on the product!** Please contact your specialist dealer in the case of greater loads or in case of any doubts. A bent or damaged walking stick must not be used under any circumstances, since it can no longer offer the necessary stability! A bent walking stick must also not be straightened because of the risk of breakage! **Note:** All serious incidents that occur in connection with the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is located.

10. Re-use (rental) / duration of use

We shall only accept liability pursuant to the Product Liability Act for the initial placing of our products on the market. We only consent to reuse if we have first inspected the product. **The service life of the walking stick is limited to 2 years.** The user shall be responsible for the continued use of the walking stick after this time.

11. Warranty and liability

You have purchased a high-quality Ossenberg product. In the event of a justified material defect in spite of the highest quality standards, Ossenberg reserves the right, as part of a 12-month warranty, to either repair the defective item or replace it with a defect-free follow-up delivery, at its own discretion. Please note that stricter legal regulations are applicable to our finished products, in particular custom-made products. We shall not be liable for any damage caused by modifications to our products. The warranty does not cover defects attributable to wear and tear or improper handling. Wear parts are excluded from the warranty. Only original replacement parts and accessories must be used! Please contact your specialist dealer for replacement parts, accessories and repairs.

12. Disposal

Please consult your specialist dealer about the regional disposal regulations.

13. Technical data

Item no. 3200

GRIP TO FLOOR

adjustable from
79 cm - 96 cm
7-fold by push button

RUBBER TIP

Item no. 970P
Inner Ø 19 mm

MAX.


Loadable up to 120 kg
body weight

WEIGHT

510 g

79 cm - 96 cm

Item no. 3201

GRIP TO FLOOR

adjustable from
80 cm - 97 cm
7-fold by push button

RUBBER TIP

Item no. 970P
Inner Ø 19 mm

MAX.


Loadable up to 120 kg
body weight

WEIGHT

400 g

80 cm - 97 cm

Item no. 3203

GRIP TO FLOOR

adjustable from
71 cm - 96 cm
11-fold by push button

RUBBER TIP

Item no. 970P
Inner Ø 19 mm

MAX.


Loadable up to 120 kg
body weight

WEIGHT

410 g

71 cm - 96 cm

14. Replacement parts

Item no. 970P Rubber tip in carbonoprik with steel inlay- inner diameter 19mm

Item no. 181P Magnetsafe Micro USB Charging Cable

15. Symbol directory

Item number

Batch number

Read the operating instructions

Attention

This product complies with the requirements of Regulation (EU) 2017/745 (MDR)

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This is a product according to the guidelines of the Medical Devices Directive. Furthermore the Ossenberg GmbH is certified according to DIN EN ISO 13485.

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Desuden er Ossenberg GmbH certificeret i